

Workspace Decontamination Log

Visby recommends daily decontamination of workspaces and immediate decontamination when a spill occurs or a suspected contamination event occurs.

To properly decontaminate the workspace, saturate a disposable paper towel with 1:10 diluted bleach (or equivalent) and wipe the surface of the workspace until it is visibly wet with the bleach solution. Let the bleach solution rest at least ten minutes on the work surface. Saturate a new paper towel with 70% ethanol to clean the bleach residue off of the work surface.

Reviewed by: _____

Date: _____

For training support please email: visby.training@visby.com. For support issues, please contact Visby Medical Customer Support at 1-833-GoVisby (1-833-468-4729) or email support@visby.com.